

PRODUCT SPECIFICATION

<b>Product name</b>	PremierPore VP
<b>Range</b>	Premier
<b>Type</b>	Sterile
<b>Marking</b>	Product name, Premier logo, code ref., lot number, size, quantity, CE mark, expiry date, latex free, sterile and single use symbols.
<b>Colour</b>	White
<b>Sterilization method</b>	Ethylene Oxide

Product code	Width (cm)	Length (cm)
5110N	6	7
5111N	10	10
5112N	15	15
5113N	10	25
5114N	10	20
5115N	10	30
5116N	10	35

<b>Primary use</b>	A dressing for open wounds. For hospital use.
--------------------	---

Packaging & Material

Product code	Rolls per inner	Rolls per case	Material
5110N	50	600	<u>Adhesive Non woven Spunlace:</u> Weight : 40g/m2 Material: 70% Rayon & 30% Polyester <u>Coating:</u> Pressure Sensitive Hot Melt Adhesive. (Polyacrylic) Weight: 30g/m2 (+ 10%) <u>Absorbent Cotton Pad:</u> Weight 130g/m2, 170/m2 Ingredient 100% Cotton <u>Release Paper:</u> Weight 80+ 10g/m2
5111N	50	600	
5112N	50	600	
5113N	50	300	
5114N	50	300	
5115N	50	300	
5116N	50	300	

<b>Quality certification</b>	ISO13485
<b>Conformity to guidelines</b>	Class 1s – Medical Device Directive 93/42/EEC
<b>Storage conditions</b>	Store in a cool dry place away from heat and direct sunlight
<b>Shelf life</b>	3 years if stored as recommended above
<b>Additional information</b>	PremierPore VP (vapour permeable) is a self adhesive absorbent island dressing that is waterproof. It is a one way stretch water repellent, non-woven dressing with a latex free SIS acrylic adhesive that provides soft, secure fixation. The wound pad is covered with a low adherent film to enable easy removal. The vapour permeable Polyurethane film coating protects the wound when showering. It also increases pad absorbency and prevents strikethrough. PremierPore is ideal for use on wounds with low or moderate levels of exudates, including post-operative wounds

*Medical Device and Reporting System:  
According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.*

<b>REVISION:</b>	6.0	<b>DATE</b>	18/02/2016	<b>AUTHOR</b>	J Bertran
<b>NOTES:</b>	Periodic review, DSN format and Shermond logo updated.				
<b>REVISION:</b>	7.0	<b>DATE</b>	29/11/2017	<b>AUTHOR</b>	F Li
<b>NOTES:</b>	Product codes are updated.				

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: <a href="http://www.shermond.com">www.shermond.com</a>	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond  
T: +44 (0) 1530 278 111  
W: [www.shermond.com](http://www.shermond.com)

Unit 6, Delta Park Industrial Estate  
Millmarsh Lane, London  
Enfield, EN3 7QJ

Shermond is a trading division of Bunzl Retail & Healthcare Supplies  
Limited Registered Office: York House, 45 Seymour Street, London W1H 7JT. Registered in England No. 6253