

**DATA SHEET: DSN084** 

## PRODUCT SPECIFICATION

Product name	Finepore Surgical Tape		
Range	Premier		
Marking	Product name, Premier logo, code ref., lot number, size, quantity, CE mark, expiry date, latex free and single use symbols.		
Colour	White		

Product code	Width (cm)	Length (m)	
5000	1.25	9.1 ± 5%	
5005	2.5	9.1 ± 5%	
5010	5.0	9.1 ± 5%	
5015	7.5	9.1 ± 5%	

Primary use	In hospital environments for holding surgical dressings and other medical devices in place on the
	patient. May also be used in other medical areas as determined appropriate by trained medical staff.

Packaging				
Product code	Rolls per	Rolls per		
	inner	case		
5000	24	240		
5005	12	120		
5010	6	60		
5015	4	40		

Material	Permeable Non-Woven Synthetic Adhesive Tape.
	The SIS Adhesive is a wholly synthetic polymer coating of ultra-low skin irritancy

Conformity to guidelines	Class 1 – Medical Device Directive 93/42/EEC		
Conformity to standards	The material meets the British Pharmacopeia specifications		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended above		
Additional information	Finepore permeable tapes are high quality surgical dressing tapes, ideal for holding tubing, surgical dressings and swaps in place.  The porous non-woven synthetic material allows moisture to escape providing excellent adhesion whilst being gentle on the skin. They are resistant to water immersion.		
	The product is non sterile and latex free.		
Medical Device and Reporting Sy			
According to the official criteria	of MDD, incidents caused by this device must be reported immediately to our offices.		

REVISON:	4.0	DATE	09/02/2016	AUTHOR	J Bertran
NOTES:	Periodic review, DSN format and Shermond logo updated.				

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond

T: +44 (0) 1530 278 111 W: <u>www.shermond.com</u> Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ