

PRODUCT SPECIFICATION
-----------------------

<b>Product name</b>	Shelf Pessary (Black)
<b>Range</b>	Universal
<b>Classification</b>	Class IIb MDD
<b>Marking</b>	Code Ref., Lot Number, Expiry date and CE Mark
<b>Material</b>	Acetylcopolymer
<b>Packaging</b>	Pouch, dispenser, carton

<b>Primary use</b>	A pessary is a small plastic or silicone medical device which is inserted into the vagina or rectum and held in place by the pelvic floor musculature
--------------------	---

Product code	Colour	Dimension mm	Quantity/outer carton
UN16179-10	Black	70	1
UN16179-11	Black	73	1
UN16179-12	Black	76	1
UN16179-13	Black	79	1
UN16179-14	Black	83	1
UN16179-15	Black	86	1
UN16179-16	Black	89	1
UN16179-17	Black	92	1
UN16179-18	Black	95	1
UN16179-19	Black	98	1
UN16179-2	Black	38	1
UN16179-20	Black	102	1
UN16179-21	Black	108	1
UN16179-22	Black	114	1
UN16179-3	Black	38	1
UN16179-4	Black	44	1
UN16179-5	Black	51	1
UN16179-6	Black	57	1
UN16179-7	Black	60	1
UN16179-8	Black	64	1
UN16179-9	Black	67	1

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: <a href="http://www.shermond.com">www.shermond.com</a>	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--

<b>Quality certification</b>	ISO 13485
<b>Storage conditions</b>	Store in a cool dry place away from heat and direct sunlight
<b>Shelf life</b>	5 Years- if stored as recommended above
<b>Additional information</b>	Shelf pessaries are single patient devices, but can be reused on the same patient after cleaning with warm soapy water. Shelf pessaries can be sterilised with Ethylene Oxide or steam at 134°C (134-138°C) for 3 mins Does not contain phthalates or latex
<i>Medical Device and Reporting System: According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.</i>	

<b>REVISION:</b>	2.0	<b>DATE</b>	03/13/2013	<b>AUTHOR</b>	H Cocking
<b>NOTES:</b>	Updated to new format and periodic review				
<b>REVISION:</b>	3.0	<b>DATE</b>	01/23/2017	<b>AUTHOR</b>	R.Zhang
<b>NOTES:</b>	Periodic review, DSN format and Shermond logo updated.				

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: <a href="http://www.shermond.com">www.shermond.com</a>	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--