

DATA SHEET: DSN324

PRODUCT SPECIFICATION

Product name	Polyfield Patient Pack
Product codes	6075AF (Small)
	6076AF (Medium)
	6077AF (Large)
Range	Premier®
Туре	Sterile
Sterilization	EO
method	



Primary use Polyfield patient packs are ideal for the use in assisting the control of infection and minimising				
	risk of cross infection during the management of the wound.			

Labelling	Product name, Premier Logo, Code Ref, lot number, Expiry Date, Single use, sterility symbol and
	CE Marking
Packing	20 packs per inner box / 80 packs per outer carton

6075AF Polyfield Patient Pack (Small)						
	PACK CONTENT					
Component Code	Description	Qty	Material			
Y20058	Blue Small Nitrile Gloves (accelerant free)	1 Pair	Nitrile			
B02213-0280	2213-0280 White Non Woven swabs 10x10cm 4ply 25g/m2 ± 10% 7		50% Polyester 50% Viscose			
Y22042	White Laminated sterile field 40x40cm	1	Laminated crepe paper			
P449428	White disposal bag 27 x 46 cm	1	PE			
Y22010	White Dressing Towel 38x43cm	1	Wood pulp fibre 43g/m2; film 15g/m2			
Y22040	Dressing Towel 38x37cm	1	Wood pulp fibre			
Y25011	White apron 69x 135cm	1	PE			
Y33020	Measuring tape 10cm	1	Paper			

6076AF Polyfield Patient Pack (Medium)						
	PACK CONTENT					
Component Code	Description	Qty	Material			
Y20059	Blue Medium Nitrile Gloves (accelerant free)	1 Pair	Nitrile			
B02213-0280	White Non Woven swabs 10x10cm 4ply 25g/m2 ± 10%	7	50% Polyester 50% Viscose			
Y22042	White Laminated sterile field 40x40cm	1	Laminated crepe paper			
P449428	White disposal bag 27 x 46 cm	1	PE			
Y22010	White Dressing Towel 38x43cm	1	Wood pulp fibre 43g/m2; film 15g/m2			
Y22040	White Dressing Towel 38x37cm	1	Wood pulp fibre			
Y25011	White apron 69x 135cm	1	PE			
Y33020	Measuring tape 10cm	1	Paper			

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond

T: +44 (0) 1530 278 111 W: <u>www.shermond.com</u> Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ



DATA SHEET: DSN324

6077AF Polyfield Patient Pack (Large)					
	PACK CONTE	NT			
Component Code	Description	Qty	Material		
Y20060	Blue Large Nitrile Gloves (accelerant free)	1 Pair	Nitrile		
B02213-0280	White Non Woven swabs 10x10cm 4ply 25g/m2 ± 10%	7	50% Polyester 50% Viscose		
Y22042	White Laminated sterile field 40x40cm	1	Laminated crepe paper		
P449428	White disposal bag 27 x 46 cm	1	PE		
Y22010	White Dressing Towel 38x43cm	1	Wood pulp fibre 43g/m2; film 15g/m2		
Y22040	White Dressing Towel 38x37cm	1	Wood pulp fibre		
Y25011	White apron 69x 135cm	1	PE		
Y33020	Measuring tape 10cm	1	Paper		

MITDI		VEC CI	SECIEL	CATI	\sim
MIIKI	LE GLO'	VES SI	CUIF	CALL	UN

Dimensions –BS EN 455-2	Small	Medium	Large	Force at break	≥ 6.0 N (during shelf life)
Palm width (mm)	85± 5	95± 5	106± 5	(BS EN 455-2)	≥ 6.0 N (after challenge
Length (mm)	240min	240min	240min		testing)

Residual powder (BS EN 455-3)	< 2 mg/glove
Residual protein	< 20 ug/g

Freedom from holes (BS EN455-1)		
Water tight test General inspection level I, AQL 1.5		

QUALITY STANDARDS

Quality certification	rality certification ISO 13485:2003				
CE Classification	Class 1s - Medical Device Directive 93/42/EEC				
Conformity to standards	The gloves confirm to BS EN455-1, 2,3 &4				
	The sterilization process is validated according to BS ISO11135				
Storage conditions	ditions Store in a cool dry place away from heat and direct sunlight				
Shelf life	helf life 3 Years if stored as recommended above				
Additional information	The product is Latex free, powder free and accelerator free. The pack remains sterile unless opened.				
Medical Device and Reporting S According to the official criteria	ystem: of MDD, incidents caused by this device must be reported immediately to our offices.				

REVISON:	6.0	DATE	17/11/2015	AUTHOR	FLi		
NOTES:	Periodic review	Periodic review, DSN format and Shermond logo updated.					
REVISON:	7.0	DATE	09/11/2017	AUTHOR	FLi		
NOTES:	Sterilization me	Sterilization method changed from Gamma to EO sterilization Swabs weight change from 30 g/ m2 to 25g/ m2					
REVISON:	8.0	DATE	18/05/2018	AUTHOR	A Gidley		
NOTES:	Updates to incl	Updates to include: addition of component codes to specifications/change of layout					

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond Unit 6, Delta Park Industrial Estate
T: +44 (0) 1530 278 111 Millmarsh Lane, London
W: www.shermond.com Enfield, EN3 7QJ