

**PRODUCT SPECIFICATION**

|                             |                     |
|-----------------------------|---------------------|
| <b>Product name</b>         | Suture Removal Pack |
| <b>Product codes</b>        | 4722                |
| <b>Range</b>                | Premier®            |
| <b>Type</b>                 | Sterile             |
| <b>Sterilization method</b> | EO                  |



|                    |  |
|--------------------|--|
| <b>Primary use</b> | For general removal of sutures.  |
| <b>Labelling</b>   | Product name, Premier Logo, Code Ref, lot number, Expiry Date, Single use, sterility symbol and CE Marking |
| <b>Packing</b>     | Each pack is contained in a paper pouch and 120bags packed into an outer carton                            |

| 4722 Suture Removal Pack |                                  |     |                                   |
|--------------------------|----------------------------------|-----|-----------------------------------|
| PACK CONTENT             |                                  |     |                                   |
| Component Code           | Description                      | Qty | Material                          |
| Y26015                   | Blue Forceps 10cm                | 1   | ABS                               |
| Y30016                   | Stitch Cutter Blade 6.5cm        | 1   | Carbon Steel                      |
| B02202-0195              | White Non Woven Swabs 5x5cm 4ply | 3   | Viscose 50% Polyester 50% 25g/ m2 |
| Y22040                   | Dressing Towel 38x37cm           | 1   | Paper                             |
| P449429                  | Orange Disposable Bag 27x46cm    | 1   | PE                                |
| Y22042                   | White Sterile Field 40x40cm      | 1   | Laminated crepe paper             |

**QUALITY STANDARDS**

|                                |  |
|--------------------------------|--|
| <b>Quality certification</b>   | ISO 13485:2003   |
| <b>CE Classification</b>       | Class Is - Medical Device Directive 93/42/EEC  |
| <b>Conformity to standards</b> | The sterilization process is validated according to BS ISO 11135                                     |
| <b>Storage conditions</b>      | Store in a cool dry place away from heat and direct sunlight   |
| <b>Shelf life</b>              | 5 Years if stored as recommended above   |
| <b>Additional information</b>  | The product is Latex free, powder free and accelerator free. The pack remains sterile unless opened. |

*Medical Device and Reporting System:  
According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.*

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

|  |  |
|--|--|
| Shermond<br>T: +44 (0) 1530 278 111<br>W: <a href="http://www.shermond.com">www.shermond.com</a> | Unit 6, Delta Park Industrial Estate<br>Millmarsh Lane, London<br>Enfield, EN3 7QJ |
|--|--|



DATA SHEET: DSN077

|                  |   |             |            |               |           |
|------------------|---|-------------|------------|---------------|-----------|
| <b>REVISION:</b> | 4.0   | <b>DATE</b> | 16/07/2015 | <b>AUTHOR</b> | F Li      |
| <b>NOTES:</b>    | Add the sterilization method and CE classification  |             |            |               |           |
| <b>REVISION:</b> | 5.0   | <b>DATE</b> | 23/07/2015 | <b>AUTHOR</b> | J Bertran |
| <b>NOTES:</b>    | Shermond Logo updated and periodic review   |             |            |               |           |
| <b>REVISION:</b> | 6.0   | <b>DATE</b> | 09/11/2017 | <b>AUTHOR</b> | F Li      |
| <b>NOTES:</b>    | Sterilization method changed from Gamma to EO sterilization<br>Swabs weight change from 30 g/ m2 to 25g/ m2 |             |            |               |           |
| <b>REVISION:</b> | 7.0   | <b>DATE</b> | 17/05/2018 | <b>AUTHOR</b> | A Gidley  |
| <b>NOTES:</b>    | Added in component codes to pack specification  |             |            |               |           |

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

|  |  |
|--|--|
| Shermond<br>T: +44 (0) 1530 278 111<br>W: <a href="http://www.shermond.com">www.shermond.com</a> | Unit 6, Delta Park Industrial Estate<br>Millmarsh Lane, London<br>Enfield, EN3 7QJ |
|--|--|

Shermond is a trading division of Bunzl Retail & Healthcare Supplies  
Limited Registered Office: York House, 45 Seymour Street, London W1H 7JT. Registered in England No. 6253