

PRODUCT SPECIFICATION

Product name	Digital Thermometers
Range	Premier
Classification	Class I MDD
Marking	Product name, Premier Logo, Size, Code Ref., Lot Number, Expiry date and CE Mark
Packaging	Pouch, dispenser, carton

Primary use	Medical thermometers are used for measuring human body temperature, with the tip of the thermometer being inserted either into the mouth under the tongue, under the armpit or into the rectum via the anus
--------------------	---

Product code	Dimension	Quantity/carton
UNE315	13x2x4cm	12

Quality certification	ISO 13485
Storage conditions	Store in a cool dry place away from heat and direct sunlight
Shelf life	5 Years- if stored as recommended above
Additional information	Does not contain phthalates or latex or mercury or glass. These digital thermometers have an easy-to-read LCD screen that is accurate to 0.1 degree centigrade and include a last temperature memory. Digital thermometers can be disinfected or covered with snug fitting probe covers. Batteries are replaceable and the devices include a low battery warning

*Medical Device and Reporting System:
According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.*

REVISION:	2.0	DATE	03/13/2013	AUTHOR	H Cocking
NOTES:	Updated to new format and periodic review				
REVISION:	3.0	DATE	01/23/2017	AUTHOR	R.Zhang
NOTES:	Periodic review, DSN format and Shermond logo updated.				

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: www.shermond.com	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--