

DATA SHEET: DSN097

Product name	Fluid protection gowns
Range	Premier
Classification	Low risk PPE
Marking	Code Ref., Description, CE Mark and our name and address
Material	Extruded polyethene
Packaging	Dispenser, carton

Primary use Staff protection, providing fluid repellence. Long sleeve, thumb loop, one size

Product code	Colour	Dimension	Qity/inner box	Qty/outer carton
5530	Blue	Length: - 120 cms (47.2 ")	20	100
		Width (body): - 96 cms (37.7 ")		
		Width (arm span) - 196 cm (77.1")		
		Thickness - 0.030 mm		
5550	Green	Length: - 120 cms (47.2 ")	20	100
		Width (body): - 96 cms (37.7 ")		
		Width (arm span) - 196 cm (77.1")		
		Thickness - 0.030 mm		
5560	Yellow	Length: - 120 cms (47.2 ")	20	100
		Width (body): - 96 cms (37.7 ")		
		Width (arm span) - 196 cm (77.1")		
		Thickness - 0.030 mm		

Storage conditions	Store in a cool dry place away from heat and direct sunlight				
Shelf life	5 Years- if stored as recommended above				
Additional information	Thumb loop fluid protection gowns are manufactured from soft durable polyethylene material of one piece construction, which is impervious to fluids, providing an extremely effective barrier against cross infection. They are suitable for staff use in areas such as isolation, delivery suites, burns units, laboratories and mortuaries. Latex Free				

Medical Device and Reporting System:

According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.

REVISON:	6.0	DATE	02/06/2017	AUTHOR	R.Zhang	
NOTES:	Periodic review, DSN format and Shermond logo updated.					

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond

T: +44 (0) 1530 278 111 W: <u>www.shermond.com</u> Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ