

PRODUCT SPECIFICATION

Product name	X-ray detectable Gauze swabs
Range	Premier
Classification	Class IIa - Medical Device
Marking	Marking on Inner Polybag consists of Product Ref.No. Our Name and Address, CE Marking, Size and Description. Marking on Outer Cases as per Inner Polybag, except for our Name and Address
Material	Pure Cotton Surgical Gauze, (Type BP13 Light) which has been hydrophilised and bleached prior to processing. The most significant feature of this material is its ability to absorb liquids rapidly with X-Ray Detectable synthetic Filament
Packaging	Polybag, carton
Sterile	Yes- ETO

Primary use	To absorb blood and other bodily fluids during surgically invasive procedures
--------------------	---

Product code	Colour	Dimension	Quantity/case	Other infor.
0104	White	10 x 7.5 cm 32 Ply	800	tied 5's
0104s	White	10 x 7.5 cm 32 Ply	160	tied 5's
0106	white	10 x 10 cm 12 Ply	70	Double wrapped
0107	White	10 x 10 cm 16 Ply	800	tied 5's
0110	White	10 x 10 cm 32 Ply	600	tied 5's
0113	White	15 x 2.5 cm 16 Ply	800	tied 5's
0116	White	15 x 10 cm 24 Ply	540	tied 5's
0122	White	36 x 11 cm 16 Ply	220	tied 5's
0124	white	35 x 22.5 cm 8ply	28	tied 5's
0130	White	22.5 x 22.5 cm 12 Ply	48	taped/tied 5's
0133	White	30 x 30 cm 12 Ply	28	taped/tied 5's
0136	White	30 x 30 cm 24 Ply	16	taped/tied 5's
0141	white	45 x 3.5 cm 12 ply	16	taped/tied 5's
0142	White	45 x 45 cm 12 Ply	12	taped/tied 5's
0144	White	53 x 22.5 cm 8ply	6	taped/tied 5's
0150	White	10 x 10 cm 2 Ply	140	tied 5's
0160	White	22.5 x 22.5 cm 12 Ply	48	tied 5's
0165	White	30 x 30 cm 12 Ply	70	tied 5's
0170	White	45 x 45 cm 12 Ply	12	tied 5's
0175	White	30 x 30 cm 18 Ply	16	tied 5's
0180	White	10 x 10 cm 32 Ply	600	tied 5's

Quality certification	ISO 13485
Conformity to standards	BP 1993
Storage conditions	Store in a cool dry place away from heat and direct sunlight

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: www.shermond.com	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--

Shelf life	5 Years- if stored as recommended above
Additional information	Latex free Each pack has 3 traceability stickers which can be removed and stuck in patient's notes.
<i>Medical Device and Reporting System: According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.</i>	

REVISION:	2.0	DATE	22/11/2012	AUTHOR	H. Cocking
NOTES:	Update to New format and Periodic Review				
REVISION:	3.0	DATE	20/09/2013	AUTHOR	H. Cocking
NOTES:	Extension of product range				
REVISION:	4.0	DATE	02/06/2017	AUTHOR	R.Zhang
NOTES:	Periodic review, DSN format and Shermond logo updated.				
REVISION:	5.0	DATE	26/07/2017	AUTHOR	J.Bertran
NOTES	Additional information included.				
REVISION:	6.0	DATE	21/08/2017	AUTHOR	R Morton
NOTES	Extension of product range				
REVISION:	7.0	DATE	29/09/2017	AUTHOR	R Morton
NOTES	Update of intended use				
REVISION:	8.0	DATE	19/10/2017	AUTHOR	J Bertran
NOTES	Amended quantities				

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Shermond T: +44 (0) 1530 278 111 W: www.shermond.com	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--